



GRACILIS®-Micro-Spoon Forceps



A fine but strong spoon forceps for safety holding of bone fragments in case of resection arthroplastics of carpus area as well as for reduction of fractures.

Through the specific formed trip with strong teeth in the front area and the lateral micro-toothing, a safety holding of smallest fragments is possible. The central depression Encloses the fragments for stabilization of handle.

A slim designed front part of this ergonomic instrument enables best visibility to surgery area.



Art.No.		Maulteil	Länge
RI.454400-12G	Spoon forceps	2.5 x 5.5 mm	120 mm

Round handle \varnothing 8mm with Gracilis-Design